

Surgical glove

# AF micro fit

**Suppressed allergens so far to protect "hands".**  
**The feeling is "Like a Second Skin". More Sensitive!**



Realized softer tactile sensitivity compared to conventional polychloroprene surgical gloves.

(Compared to our company products)

AF microfit is a new synthetic rubber surgical gloves that do not contain casein, powder and natural rubber latex as well as accelerators.

\* Feature of polychloroprene surgical gloves is more stable durability compared to conventional surgical gloves (made of natural rubber latex and polyisoprene synthetic rubber).

## Accelerator Free

**"Accelerators"** is said to occupy 80%\* of allergic contact dermatitis caused by gloves. Also be careful with synthetic rubber gloves, as some products contain accelerators.

\*Reference:

Allergic and irritant reactions to rubber gloves in medical health services.  
Journal of the American Academy of Dermatology Volume25,Number 5,  
Part1 November 1991

## FEATURE

### ● Aim for a better fit

For medical professionals who require delicate hand sensation, achieve a thin design and ensure excellent wear ability and grip by Ansell's own manufacturing method.

### ● Eliminate as much risk as possible due to surgical glove

- Accelerator Free
- Powder Free
- Latex Free (made of polychloroprene synthetic rubber)

## SPECIFICATION

Color	Cream
Thinness (Fingertip)	0.17 mm
Thinness (Palm)	0.16 mm
Thinness (Cuff)	0.15 mm
Shape	Curved finger
Internal coating	Smooth

\*In-house measured results

Product Number	Size	Package
805531710	5.5	50 Pairs × 4 Boxes = 200 Pairs
806031710	6	
806531710	6.5	
807031710	7	
807531710	7.5	
808031710	8	
808531710	8.5	

GAMMEX is a trademark of ANSELL LTD.