

Examination glove

# Nitrile Glove

## Non-natural rubber Examination and Treatment Use Rubber Glove



### FEATURE

- **Latex-free**

Does not contain natural rubber latex causing Type I allergies.

- **Powder-free**

Reduce the risk of complications and dermatitis due to free of powder.

- **Excellent fingertip feeling & fit**

It is a thin type emphasizing fingertip feeling and is stretchable.

Therefore, it fits the hand and can be used for a wide range of procedures.

- **Safety - Protection from drug exposure**

Anticancer drug penetration test (ASTM D6978-05) has been conducted.

- **Embossing to prevent tape from sticking**

It is also suitable for blood collection by embossing on the fingertips.

- **Cuff to prevent flip-up**

Beaded cuff prevents from flip-up.

- **Clean color**

Adopt white as clean image color.

## **!! Caution for use !!**

- This product is a medical disposable product and should be used only once.
- Stop to use if symptoms such as itching, rash, or something appear.  
(Depending on the constitution, the above symptoms may occur.)
- Do not use the product if it is found to be damaged, deformed, or has pinholes or other abnormalities.
- Be sure to check before use as some medicine and pharmaceuticals may swell or permeate.
- It may tear or puncture if you touch a sharp tool such as a knife. Please check before use.
- If this product tears during use, replace it with a new one.
- Dispose in a safe way with consideration for preventing infection after use.
- Please read the attached document before use.
- ASTM is an international industrial standard developed and published by ASTM International (formerly known as the American Society for Testing and Materials (ASTM)).  
Please contact us for details of test results on this product.

## **SPECIFICATION**

Product Number	Size	Package
809059540	SS	200 Pcs × 10 Boxes = 2,000 Pcs
809069540	S	
809079540	M	
809089540	L	