

Surgical glove

Micro-Thin Powder Free UG

Safe and secure utility glove for all purposes



FEATURE

- A polymer coating of polyurethane is applied to the inner surface for easy wearing for intraoperative exchange, double fitting, or with wet hands.
- Adopt brown color to prevent glare caused by reflections under surgical lights and microsurgery.
- Suitable for microsurgery due to very thin and excellent tactile sensitivity.
- To minimize the risk of sensitization of latex allergy and chemical allergy, eliminate latex proteins and chemicals as much as possible by our own special cleaning method.
- Reduce pressure during double fitting due to original shape.
- Prevent cuff roll-down by using Sure Fit technology.

!! Caution for use and storage!!

- This product is made of natural rubber and should not be used by medical professionals or patients who have experienced anaphylactic-like symptoms in the past. Natural rubber may rarely cause allergic symptoms such as itching, redness, hives, swelling, fever, difficulty in breathing, symptoms like asthma, hypotension, and allergic symptoms such as shock.
Allergic symptoms can be Type I (immediacy: due to proteins included in natural rubber latex) or Type IV (deferred: mainly due to chemicals of additives such as accelerators).
- Be careful regarding the use of this product by those who have frequent exposure to natural rubber, such as patients with spina bifida as they are considered a high-risk group who may occur allergic symptoms of natural rubber.
- This product is for single use only and should not be reused.
- This product should always be used under the supervision of a doctor.
- Be careful not to get wet, and store this product away from high temperatures, high humidity, direct sunlight, or UV rays such as fluorescent lamps.
- Do not use if the packaging is damaged or stained, and if this product is found to be damaged, deformed or other abnormalities such as pinholes.

SPECIFICATION

Product Number	Size	Package
805531610	5.5	50 Pairs × 4 Boxes = 200 Pairs
806031610	6	
806531610	6.5	
807031610	7	
807531610	7.5	
808031610	8	
808531610	8.5	

GAMMEX is a trademark of ANSELL LTD.



JMS CO.,LTD.
12-17 KAKO-MACHI, NAKA-KU,
HIROSHIMA 730-8652, JAPAN
TEL: +81-(0)82-243-5505
<https://www.jms.cc/>

2024.03.000H003-JMS