




INFUSION SET



For safe infusion therapy

Are you using a lock type of infusion set now?

The luer type has a weak connection and is easily detached when an external force is applied.

Compared to that, the lock type cannot be easily detached due to the screwed type. For safe infusion therapy, use a lock type of infusion set.

Infusion Set (20 drops/mL)

Product Code	Name	Tip Shape	Quantity
014502043	INFUSION SET 20 JY220 ADAPTA DEHPF	Luer Fit	400 pcs
014502042	INFUSION SET 20 JY220 LOCK DEHPF	Luer Lock	400 pcs

Infusion Set (60 drops/mL)

Product Code	Name	Tip Shape	Quantity
018802042	INFUSION SET PDIA 20JY220 ADAPTA DEHPF	Luer Fit	400 pcs
018802040	INFUSION SET PEDIA JY220 LOCK DEHPF	Luer Lock	400 pcs



JMS Co., Ltd.

12-17 Kako-machi, Naka-ku, Hiroshima 730-8652, Japan

TEL: +81-(0)82-243-5505 Fax: +81-(0)82-243-1976

<http://www.jms.cc>

2024.03.000H015-JMS