









INJECTION NEEDLE



Features :

- Super thin-wall needle allows for one size larger flow rate even with the same outer diameter.
- Individual packaging is made of plastic, eliminating the need for separate disposal.
- Separate packaging eliminates the need to separate the product for use.
- The needle base is translucent, making it easy to check the blood flush bag during intravenous injection.

Injection Needle Line-up

18G pink (1.2mm)	 1½" (38mm)  1½" (38mm)	23G deep blue (0.6mm)	 1¼" (32mm)  1¼" (32mm)  1" (25mm)
19G cream (1.1mm)	 1½" (38mm)  1½" (38mm)	24G medium purple (0.55mm)	 1¼" (32mm)  1" (25mm)
20G yellow (0.9mm)	 1½" (38mm)  1½" (38mm)	25G orange (0.5mm)	 1" (25mm)
21G deep green (0.8mm)	 1½" (38mm)  1½" (38mm)	26G brown (0.45mm)	 ½" (13mm)  ½" (13mm)
22G black (0.7mm)	 1½" (38mm)  1½" (38mm)  1¼" (32mm)  1¼" (32mm)	27G medium grey (0.4mm)	 ¾" (19mm)  ½" (13mm)

Product Code	Product Name	Tip shape	Quantity
501864102	INJECTION NEEDLE 18GX11/2 RB	Regular bevel	4000 pcs
501864202	INJECTION NEEDLE 18GX11/2 SB	Short bevel	4000 pcs
501934101	INJECTION NEEDLE 19GX3/4 RB	Regular bevel	4000 pcs
501964102	INJECTION NEEDLE 19GX11/2 RB	Regular bevel	4000 pcs
502024101	INJECTION NEEDLE 20GX5/8 RB	Regular bevel	4000 pcs
502064102	INJECTION NEEDLE 20GX11/2 RB	Regular bevel	4000 pcs
502064202	INJECTION NEEDLE 20GX11/2 SB	Short bevel	4000 pcs
502124102	INJECTION NEEDLE 21GX5/8 RB	Regular bevel	4000 pcs
502144102	INJECTION NEEDLE 21GX1 RB	Regular bevel	4000 pcs
502164102	INJECTION NEEDLE 21GX11/2 RB	Regular bevel	4000 pcs
502164202	INJECTION NEEDLE 21GX11/2 SB	Short bevel	4000 pcs
502244102	INJECTION NEEDLE 22GX1 RB	Regular bevel	4000 pcs
502254102	INJECTION NEEDLE 22GX11/4 RB	Regular bevel	4000 pcs
502254202	INJECTION NEEDLE 22GX11/4 SB	Short bevel	4000 pcs
502264102	INJECTION NEEDLE 22GX11/2 RB	Regular bevel	4000 pcs
502264202	INJECTION NEEDLE 22GX11/2 SB	Short bevel	4000 pcs
502344102	INJECTION NEEDLE 23GX1 RB	Regular bevel	4000 pcs
502344202	INJECTION NEEDLE 23GX1 SB	Short bevel	4000 pcs
502354102	INJECTION NEEDLE 23GX11/4 RB	Regular bevel	4000 pcs
502354202	INJECTION NEEDLE 23GX11/4 SB	Short bevel	4000 pcs
502444102	INJECTION NEEDLE 24GX1 RB	Regular bevel	4000 pcs
502454102	INJECTION NEEDLE 24GX11/4 RB	Regular bevel	4000 pcs
502524102	INJECTION NEEDLE 25GX5/8 RB	Regular bevel	4000 pcs
502544102	INJECTION NEEDLE 25GX1 RB	Regular bevel	4000 pcs
502614102	INJECTION NEEDLE 26GX1/2 RB	Regular bevel	4000 pcs
502614202	INJECTION NEEDLE 26GX1/2 SB	Short bevel	4000 pcs
502714102	INJECTION NEEDLE 27GX1/2 RB	Regular bevel	4000 pcs
502714202	INJECTION NEEDLE 27GX1/2 SB	Short bevel	4000 pcs
502734202	INJECTION NEEDLE 27GX3/4 RB	Regular bevel	4000 pcs



JMS Co., Ltd.

12-17 Kako-machi, Naka-ku, Hiroshima 730-8652, Japan

Phone: +81-(0)82-243-5505 Fax: +81-(0)82-243-1976

<http://www.jms.cc>

2024.03.000H021-JMS

2025.03.000H021-JMS