

Transfusion Set

FEATURE

As this product uses a large sized filter, it reduces the loss of blood components and maintains the required transfusion rate

- ◆ Filter mesh size: 100 mesh (approx. 180μm).
- ◆ The clamp is red to easily identify as a transfusion set.

SPECIFICATION

Standard (20 drops/mL)

Product Number	Tip	Spike	Package	Remarks
041140301	Lock	Middle tube type	25 Pcs × 10 Boxes = 250 Pcs	
041040200	Lock	Middle tube type	25 Pcs × 10 Boxes = 250 Pcs	Three way stopcock × 1
041040100	Lock	Middle tube type	25 Pcs × 10 Boxes = 250 Pcs	Three way stopcock × 2
043540400	Lock	Direct type	25 Pcs × 10 Boxes = 250 Pcs	



Type Y (20 drops/mL)

Suitable for massive blood transfusions.

Product Number	Tip	Package	Remarks
049040101	Lock	25 Pcs × 10 Boxes = 250 Pcs	
049040200	Lock	25 Pcs × 10 Boxes = 250 Pcs	Three way stopcock × 1

