

For Hemodialysis

Blood Tubing Set

Support a safe and more secure medical environment

FEATURE

- ◆ JMS Needle-less Access Port: Planecta™
Reduce the risk of infection and needlestick accidents
- ◆ JMS spiral chamber
Improve blood flow conditions in the air trap chamber

Cannula connector cap

- Recirculation cap
 - AL/BL/VL Cap
- ① Securely cover the cannula connector part
 - ② AL/VL can be shorted without removing the cap
 - ③ Secure more cleanly the cannula connector part

Tightly secure clamp

- Large clamp that tightly closes the main clamp
- Compact clamp f that tightly closes the subline

Cannula connector

Available for free lock and slide lock

Air trap chamber

Select the size and type according to the purpose of use

Luer lock connector

Enable to be connected to Luer-lock syringe

Planecta™

Planecta™ does not use a puncture needle

- ① Infection prevention
- ② No dead space

